

Very Long Path Length Ion Mobility Separations using Structures Lossless Ion Manipulations

Richard D. Smith, Ian K. Webb, Ahmed Hamid, Sandilya Garimella, Yehia M. Ibrahim,
Aleksy V. Tolmachev, Spencer A. Proust, Gordon A. Anderson and Erin S. Baker


Pacific Northwest National Laboratory

Some approaches for potentially achieving very high resolution ion mobility separations


- Separation over long path in dense gas/fluid; challenges: large voltage drop required, difficulty in combining with MS
- Separation in very low temperature fluid; challenges: design complexity, maintaining uniform temperature over useful volume
- Extended residence time Differential Mobility separations (e.g. FAIMS); challenges: large ion losses, slow scan speed
- Separation in flowing (or expanding gas); challenges: maintaining highly stable gas flow, limited separation space and dynamic range
- Separation in cyclic path devices; challenge: limited separation space
- **Separation over very long path lengths; challenges: space and cost needed for long path device, large voltage drop required**

Structures for Lossless Ion Manipulations (SLIM)

One of many SLIM electrode arrangements.....



See Wednesday Poster by Tolmachev et al., WP 485




Garimella et al., JASMS, 25, 1890 (2014)

A traveling wave SLIM IMS variation


- Existing data shows IMS performance can approach drift tube designs, and theory suggests Resolution \propto (drift path length)^{0.5} with good separations achievable over significant mobility range*
- Attraction: voltages applied independent of drift path length
- Question: can a SLIM implementation provide good resolution?
- 2nd question: can the SLIM implementation be simplified?
- 3rd question: can the design be made compact?
- 4th question:

* Shvartsburg and Smith, Anal. Chem., 80, 9689 (2008)

Initial traveling wave IMS SLIM electrode design evaluated




#RF strips, #rows DC electrodes (12,11)



5 mm PCB gap

Elevated DC to 4 of every 8 electrodes in each of 11 rows incremented in single electrode steps




Sequence = 11110000


Ion transmission efficiency by trajectory simulations

Board spacing: 5 mm gap

Guard 15 V, TW amplitude 30 V and speed 84 m/s




Experimental arrangement



Direct measurement of ion current in SLIM IMS module


Board spacing: 5 mm

Guard 15V, TW amplitude 30 V and speed 84 m/s



See Monday Poster by Hamid et al.; MP 137


TW speed effect on IMS resolution




Sequence = 11110000

Resolving power ~ 35


10.5 m/s



84 m/s




210 m/s




See Monday Poster by Hamid et al.; MP 137

Good agreement between experiment and simulation




Comparison of 30 cm TW SLIM with conventional TW IMS

**TW SLIM
(30 cm)**



- P = 4.0 Torr N₂
- TW speed = 89 m/s
- TW amp. = 30 V

**Synapt G2-TWIMS*
(25 cm)**




- P = 2.6 Torr N₂
- TW speed = 650 m/s
- TW amp. = 40 V


Liulin Deng

*Li, H.; Bendiak, B.; Siems, W. F.; Gang, D. R.; Hill Jr, H. H. *Anal. Chem.*, 85, 2760 (2013)


A few of the TW sequences examined in SLIM modules



Effect of TW sequence on IMS resolution for 6,5 arrangement with 5 mm board gap



TW amplitude effect on IMS resolution



Sequence = 11110000

Some other SLIM electrode arrangements explored

#RF,#DC

5,6


3,4

4,3

2,3



3,2

3,2




IMS performance for a 6,5 electrode arrangement

#RF strips, #rows DC electrodes (6,5)



30.5 cm

Effect of SLIM board gap on R for 6,5 linear path arrangement



5 mm gap

Resolution = 9.8



3 mm gap


Resolution = 11.2



2 mm gap

Resolution = 9.3

Effect of SLIM board gap on R for 6,5 multi-turn path arrangement



5 mm gap

Resolution = 10.8



3 mm gap



Resolution = 11.6



2 mm gap

Resolution = 9.8


Resolution for straight vs. multi-turn path arrangements



See Monday Poster by Hamid et al.; MP 137

Simple electrode arrangement for turning corners


Ion trajectories calculated for 3 mm gap, TW amplitude = 30 V, speed = 84 m/s



The long or winding roads towards higher resolution




28 M path SLIM IMS module incorporating ~60,000 discrete electrodes




Initial evaluation of long path length SLIM IMS-MS

- ~65% efficiency for ion transmission through long path module
- Guard electrode spacing between tracks too narrow; ions jump paths at turns* under conditions needed for separation
- Refined design developed





Take 3 long path (14 M) SLIM design

TW Electrode dimensions: 0.43 mm width x 0.91 length; 0.13 electrode gaps



Initial results....much optimization to come



A good separation is often not enough


Fast switching of selected ion packets to different path previously demonstrated (Webb et al., Anal. Chem., **86**, 9169; 2014)

Trapping and accumulation of selected packets from multiple separations (see Tuesday Poster by Tsung-Chi Chen et al.; TP 074)

Our plan: Capture entire separations in a large array of traps, with intensity dependent number of accumulation steps for each trap (i.e., AGC) to increase dynamic range


Goal: Enable 'read out' of the entire separation, or **any** selected fractions, at **any** speed, to **any** MS platform type

Screen shot of ion trajectory simulation; ions of three different mobilities first separated and shown fractionated into three different traps, stored, and selected (green) for transfer to MS



Acknowledgements

PNNL SLIM Team



Xinyu Zhang
Brian LaMarche
Randy Norheim
Marques Harrer
Keqi Tang
Jonathan Cox

PNNL Laboratory Directed R&D

National Institute of General Medical Sciences

Biomedical Technology Research Resource at PNNL

DOE Office of Biological and Environmental Research