

# A carrier-assisted targeted mass spectrometry approach for proteomics analysis of single cells

Tujin Shi<sup>1</sup>, Matthew J. Gaffrey<sup>1</sup>, Thomas L. Fillmore<sup>2</sup>, Song Nie<sup>1</sup>, Carrie D. Nicora<sup>1</sup>, H. Steven Wiley<sup>2</sup>, Karin D. Rodland<sup>1</sup>, Tao Liu<sup>1</sup>, Richard D. Smith<sup>1</sup>, and Wei-Jun Qian<sup>1</sup>

<sup>1</sup>Integrative Omics Group, Biological Sciences Division, Pacific Northwest National Laboratory, Richland, WA

<sup>2</sup>Environmental Molecular Sciences Laboratory, Pacific Northwest National Laboratory, Richland, WA



Pacific Northwest  
NATIONAL LABORATORY

Proudly Operated by Battelle Since 1965

## Overview


- A simple, convenient targeted MS approach, cPRISM-SRM, was developed for enabling highly sensitive quantification of proteins in single cells.
- Proteins with  $\geq 25,000$  copies/cell in 10 mammalian cells can be reliably detected by cPRISM-SRM with minimal sample loss.
- cPRISM-SRM opens a new avenue in biomedical research involving single cells or extremely small amounts of precious clinical specimens.

## Introduction


- Antibody-based flow cytometry and mass cytometry are predominant technologies for targeted proteomics analysis of single cells.
- Cytometry methods share common shortcomings with other antibody-based methods (e.g., low-multiplex, the need of high-quality antibodies, and unavailability of antibodies for new proteins), as well as lack quantitation accuracy to provide accurate protein concentrations.<sup>1,2</sup>
- MS-based targeted proteomics has emerged as a promising alternative for antibody-free, precise, high-multiplex quantification of target proteins.<sup>3</sup>
- There is an unmet challenge for MS-based proteomics analysis of single cells because of ineffective sampling and insufficient MS sensitivity.<sup>4</sup>

## Methods

### cPRISM-SRM



### PRISM-SRM<sup>5,6</sup>




**Figure 1.** cPRISM-SRM workflow for proteomics analysis of single cells. Exogenous BSA protein was used as a carrier for lossless processing of single cells. A highly sensitive targeted MS platform PRISM-SRM was then used for absolute quantification of target proteins.

### Two types of a small number of cells


- Mimic cells: an equivalent small number of cells from large HMEC digests spiked into  $\sim 25 \mu\text{g}$  of BSA digests.
- Intact cells: a small number of intact HMEC cells (from serial dilution) collected into a tube containing  $50 \mu\text{g}$  of BSA followed by normal sample processing.

## Results

### Mimic HMEC cells




### Calibration curves



**Figure 3.** Calibration curves for quantification of K/NRAS (SYGIPFIETSAK) and SHC1 (EASLVCEAVPGAK) at 1-1000 cells.

### Intact HMEC cells (10 vs 100 cells)



**Figure 4.** XICs of transitions monitored by cPRISM-SRM at 10 and 100 intact HMEC cells isolated from serial dilution.

## Conclusions

- We developed a new targeted MS approach, cPRISM-SRM, that couples carrier-assisted sample preparation (to address ineffective sampling) to a highly sensitive targeted MS platform (to address MS sensitivity issue) for enabling proteomics analysis of single cells.
- cPRISM-SRM was demonstrated for reliable quantification of EGFR pathway proteins at  $\geq 25,000$  copies/human cell in single cells (e.g., 10 HMECs).
- We envision that this new targeted MS capability will have broad utility in biomedical research for absolute protein quantification in small-size samples that are not readily accessible by current proteomics techniques (e.g., small subpopulations of cancer cells and circulating tumor cells).

## Acknowledgements

Portions of the research were supported by U24-CA-160019 (to RDS) and U24-CA-210955 (to TL) from the National Cancer Institute Clinical Proteomic Tumor Analysis Consortium (CPTAC), and P41GM103493 (to RDS). Samples were analyzed using capabilities developed under the support of NIH National Institute of General Medical Sciences (8 P41 GM103493-10). This project was supported by the U.S. Department of Energy (DOE) Office of Biological and Environmental Research (OBER) Pan-omics program at Pacific Northwest National Laboratory (PNNL) and performed in the Environmental Molecular Sciences Laboratory, a DOE OBER national scientific user facility on the PNNL campus. PNNL is a multiprogram national laboratory operated by Battelle for the DOE under contract DE-AC05-76RL01830.

## References

- Spitzer MH, Nolan GP. *Cell* 165:780-791 (2016).
- Finck R, Simonds EF, Jager A, Krishnaswamy S, Sachs K, Fantl W, Pe'er D, Nolan GP, Bendall SC. *Cytometry A* 83:483-494 (2013).
- Shi T, Su D, Liu T, Tang K, Camp DG 2nd, Qian WJ, Smith RD. *Proteomics* 12:1074-1092 (2012).
- Wang D, Bodovitz S. *Trends Biotechnol* 28:281-290 (2010).
- Shi T, Fillmore TL, Sun X, Zhao R, Schepmoes AA, Hossain M, Xie F, Wu S, Kim JS, Jones N, Moore RJ, Pasa-Tolić L, Kagan J, Rodland KD, Liu T, Tang K, Camp DG 2nd, Smith RD, Qian WJ. *Proc Natl Acad Sci USA* 109:15395-15400 (2012).
- Shi T, Sun X, Gao Y, Fillmore TL, Schepmoes AA, Zhao R, He J, Moore RJ, Kagan J, Rodland KD, Liu T, Liu AY, Smith RD, Tang K, Camp DG 2nd, Qian WJ. *J Proteome Res* 12:3353-3361 (2013).
- Shi T, Niepel M, McDermott JE, Gao Y, Nicora CD, Chrisler WB, Markillie LM, Petyuk VA, Smith RD, Rodland KD, Sorger PK, Qian WJ, Wiley HS. *Sci Signal* 9:rs6 (2016).

**CONTACT:** Tujin Shi, Ph.D.  
Biological Sciences Division  
Pacific Northwest National Laboratory  
E-mail: [tujin.shi@pnnl.gov](mailto:tujin.shi@pnnl.gov)

**Career Opportunities:** For potential openings in Integrative Omics Group at PNNL please visit <http://omics.pnl.gov/careers>

[www.omics.pnl.gov](http://www.omics.pnl.gov)

@omicsPNNL

